

# Snap Medical kit comparison





**EpinephrineSnap®-V**



**EpinephrineSnap®-EMS**



**EpinephrineSnap®**

	EpinephrineSnap®-V	EpinephrineSnap®-EMS	EpinephrineSnap®
<b>NDC number</b>	71923-200-20	71923-100-20	71923-300-13
<b>Epinephrine Presentation</b>	1 vial epinephrine 1 mg/ mL 1:1000	1 vial epinephrine 1 mg/ mL 1:1000	1 vial epinephrine 1 mg/ mL 1:1000
<b>Number of Doses</b>	up to <b>3</b>	up to <b>2</b>	up to <b>3</b>
<b>Indication</b>	Patients weighing equal to or greater than 15 kg.	Patients weighing equal to or greater than 15 kg.	Pediatric and adult patients
<b>Needles</b>	<b>(3)</b> 23-gauge 1-inch Safety Needles	<b>(2)</b> 23-gauge 1-inch Safety Needles	<b>(3)</b> 23-gauge 1-inch Safety Needles
<b>Syringes</b>	<b>(3)</b> syringes 	<b>(2)</b> syringes	<b>(3)</b> syringes 
	Our proprietary <b>The Epinephrine Syringe™</b> is an FDA approved 1 mL syringe that has markings for common epinephrine dosing of 0.15 mg or 0.15 mL, 0.3 mg or 0.3 mL, and 0.5 mg or 0.5 mL. Always follow your prescribed order when determining epinephrine dosing and if this product will meet your protocol for treating anaphylaxis.		The 1 mL <b>Non-Residue Gasket For Minimal Dead Space Syringe</b> provides for more precise dosing, for even the littlest of patients while reducing drug wasted.
<b>FDA-approved Repackaging Facility</b>	Yes	Yes	Yes
<b>FDA Status</b>	Registered	Registered	Registered
<b>Pharmacy Board Licensed Manufacturer</b>	Yes	Yes	Yes

**Value** The team at Snap Medical Industries spent two years researching and vetting ideas for our epinephrine products based on clinical experience and use of epinephrine in anaphylaxis. We are a licensed manufacturer of epinephrine convenience kits with the boards of pharmacy nationwide. Our products are registered with the FDA and we retain FDA experts to lead our Quality Management System. These steps provide assurances that EpinephrineSnap® kits are manufactured under the highest standards available worldwide.

**Choice** EpinephrineSnap® kits are an alternative for healthcare providers, to high cost and largely unavailable epinephrine auto-injectors.

**Indication for Use** EpinephrineSnap®-V and EpinephrineSnap®-EMS can be used to treat anaphylaxis in people weighing equal to or greater than 15 kg or 33 lbs. The size of an average 2 year old according to CDC published growth charts.

EpinephrineSnap® can be used for pediatric and adult patients.

**Product and Packaging** EpinephrineSnap® kits include FDA approved epinephrine (vial) USP 1:1000 1 mg/ mL with 1 mL total volume. Our packaging is sturdy and provides no chance of light to penetrate, protecting the effectiveness of the epinephrine.

**Safety Needles** We provide 23 gauge 1 inch safety needles, assuring an easy, quick fill from the vial. The needle size is appropriate for a pediatric, teen or adult intra-muscular injection into the vastus lateralis muscle, the preferred injection site.

**Syringes** Both types of syringes have a luer lock tip. This is a critical factor as the locking action of the needle to the syringe allows the health provider to confidently administer lifesaving epinephrine in the field, even thru clothing, and never worry about the needle "slipping".

**Instruction for Use** Our kits are labeled clearly and provide the necessary information to treat anaphylaxis. Our Instructions for Use have gone thru a vigorous quality review and serve as both instructions and a reference for providers.

Learn more at [snapmedicalindustries.com](http://snapmedicalindustries.com)

EpinephrineSnap®-V and EpinephrineSnap®-EMS and The Epinephrine Syringe™ are Patent Pending.